



---

# **Block 1 Training Tools Installation Guide**

Document Publish Date: 15 June 2004

Document Version 1



Change History			
Date	Name	Change Description	Source:
5/11/2004	Steve Patterson	Created	
9/2/2004	Steve Patterson	Updated the Instructional Video icon	

## Table of Contents

Training Tools Instructional Video .....	5
Medical Web Tutorials .....	6



---


# Training Tools Instructional Video

## Hardware/Software Requirements

- Operating Systems:
  - Windows 2000 Service Pack 2
  - Windows NT 4.0 Service Pack 6
  - Windows XP Professional SP 1
- Process: Pentium II 300Mhz (Minimum)
- Minimum of 128 MB of RAM, and at least 500 MB free of disk space
- Computer Screen Resolution: 800 x 600 (Minimum)
- Color Depth: Thousands Colors – 16-bit (Minimum – Millions of Colors – 24 or 32-bit (Recommended)
- CD-ROM: 8x CD-ROM Drive (Minimum)
- Video: MPEG-1 capability
- Audio: Sound capability
- Web Browser supported:
  - Internet Explorer (IE) 5.01, 5.5 and 6.0
  - Netscape 4.7.x

## Installations Instructions

The CD has an autorun feature. When you insert Training Tools Instructional Video disc #1 CD into your CD-ROM drive, the InstallShield Wizard appears. Accept the InstallShield Wizard application defaults and the install process will load all necessary files including the following desktop icon.



If the autorun feature fails to launch, navigate to your CD drive and double-click **Setup.exe** to begin the InstallShield Wizard.

**Note:** If a previous version of the Training Tools Offline Medical Videos is installed, it is recommended the previous version be uninstalled.

---

# Medical Web Tutorials

## Hardware/Software Requirements

- Computer Processor: Pentium 166Mhz
- Memory: 64 MB (Minimum)
- Computer Screen Resolution: 800 x 600
- Operating Systems:
  - Microsoft Windows NT 4.0 Service Pack 6
  - Microsoft Windows 98
  - Microsoft Windows ME
  - Microsoft Windows 2000 Service Pack 2
- Web Browsers supported:
  - Internet Explorer 5.01, 5.5 and 6.0
  - Netscape 4.7.x
- Java Virtual Machine (JVM):
  - JVM 5.0.3085 and higher for IE
  - JVM 1.1.3 and higher for Netscape
  - Sun J2RE 1.4.2\_01
- Microsoft Media Player 6.4, 7, and 9

## Installation Instructions

The CD has an autorun feature. When you insert Medical Web Tutorials (Training Tools Offline Website Block 1) CD into your CD-ROM drive, the InstallShield Wizard appears. Accept the InstallShield Wizard application defaults and the install process will load all necessary files including the following desktop icon be installed.



If the autorun feature fails to launch, navigate to your CD-ROM drive and double-click **Setup.exe** to begin the InstallShield Wizard.

**Note:** If a previous version of the Training Tools offline Website Block 1 is installed, it is recommended the previous version be uninstalled.